Colorado (P)

Plasma, Planetary Surfaces, and Cosmic Dust Experiments at the University of Colorado



Tobin Munsat University of Colorado

Bringing Space Down to Earth Workshop, UCLA April 10, 2017

(SEPVI

IMPAG

impact.colorado.edu





Experimental setup and parameters							
				lunar casa			
		Parameter	Laboratory	(Strong B region)			
Plasma flow	Insulating plate	lon species	N ₂ +	H+			
		Ion flow energy E_b (eV)	100 - 800	1000			
		Electron Temperature T_e (eV)	0.5 (cold), 10 (hot)	10			
		Ion Temperature T _i (eV)	14	10			
		Ion Mach number M	11	9			
		Electron gyro ratio (r _e / L)	< 1 (0.3 cm / 2 cm)	<< 1 (0.35 km / 30 km)			
		Ion gyro ratio (r _i / L)	>> 1 (250 – 720 cm / 2 cm)	> 1 (150 km / 30 km)			
	B field	Electron Debye ratio (λ_{De} / L)	< 1 (0.2 cm / 2 cm)	<< 1 (0.01 km / 30 km)			
		lon Debye ratio (λ_{Di} / L)	< 1 to > 1 (1.4 – 6.5 cm / 2 cm)	<< 1 (0.1 km / 30 km)			
B field strength ~580G maximum		* The magnetic field strength 30 nT at 30 km altitude is used for the lunar case with strong magnetic anomalies [<i>Hood et al.</i> , 2001].					
impact.colorado.edu SERM IMPACT Colorado @ 1							

Dust is Ubiquitous in the Universe						
	Star Formation	Planet Formation	Solar System Evolution			
impact.color	ado.edu	NAME SERVI	IMPACT Colorad	O 😥 26		

Ablation of Micrometeoroids in Earth's Atmosphere

6 Mass influx / 10⁹ g yr⁻¹ mass decade⁻¹ Every day, billions of micrometeoroids IDPs Huge impactors ablate in Earth's upper atmosphere 5 Radars The ablated materials affect a wide • variety of atmospheric phenomena (metal layers, noctilucent clouds, etc.) 3 The micrometeoroids can inform us Visual meteors Satellite 2 detectors about the dust environment of the Earth 1 Earth's atmosphere can be used a dust Meteorites detector 10-10 10-15 1010 10-5 10⁰ 105 1015 Interpretation of radar signals depends • Mass of particle / g on knowing ionization coefficient (β) [J. Plane, Chem. Soc. Rev. 2012] Colorado (P) impact.colorado.edu 31

Gas Target/Chamber for Ablation Studies

- Differentially-pumped target capable of up to 0.5 Torr
- Complete / partial ablation of within 50 cm chamber
- Series of amplified electrodes with bias
- Measure Charge Production
- Series of viewing ports, fibercoupled optics, PMT's:
- Measure Light Production

impact.colorado.edu

Cryogenic Target for Studies of Ice Impacts

- Cryogenic system based on LN₂ flow
- · Embedded heaters control temperature baseline
- Target temperature ≥ 90K
- Two Separate Ice Targets:
 - Frozen Liquid Target
 - Deposited Thin Films
- UHV Chamber (oil-free vacuum systems)
- Ultimate Capabilities:
 - Crystalline Ice
 - Ice / Dust mixtures

impact.colorado.edu

38

Colorado 😰

Analysis Tools for Ice Growth and Impact Products

- Fabry-Perot interferometer for measuring ice thickness during deposition
- Linear time-of-flight mass-spectrometer to assess impact products
- Target substrate biased up to ±3.5 kV
- Grounded grid & target, 1 m flight tube

MCF

Shock Chemistry Studies with Mass-Spec Data Co-added CDA spectra of salt-rich ice grains • Testing hypotheses of synthesis of (Na2CO3)Na⁺ (NaCl)Na Log Amplitude complex organic materials based on micrometeoroid energy + material Comparing to previous laboratory ٠ Na⁺ (H20)Na⁺ (NaOH)Na⁺ . (NaOH)2Na⁺ . NaOH)₃Na studies using lasers and/or other radiation forms Log Amplitude Interpretation of mass spectra from • SUDA instrument on Europa mission Salt water calibration spectrum Calibration of SUDA dust detector Time of Flight • [F. Postberg, Nature 2009] Colorado (P) impact.colorado.edu 42 NASA